

## Molecular Analysis of Nanoscale Defects and Residue

NCCAVS TFUG: Advances in Metrology July 28, 2021 sung@molecularvista.com

## Molecular Vista Inc (Introduction to the Company)





## Current Status of Analytical Techniques

	Raman	FTIR	TOF-SIMS	XPS	TXRF	SEM/EDS	TEM	Auger
Species Detected	M.I.	M.I.	M.I.	M.I.	E.I.	E.I.	E.I.	E.I.
Imaging/Ma pping	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Lateral Resolution	>0.5 μm	> 10 µm	>0.2 μm	10 μm – 2 mm	~ 10 mm	1 nm* 0.5μm EDS	0.2 nm	> 10 nm
Depth Probed	> 500 nm	1 µm	1 nm	10 nm	10 nm	1 µm	~ 100 nm	10 nm

\* Imaging

M.I. Molecular information

on E.I. Elemental information



## Infrared Photo-induced Force Microscopy (IR PiFM)



Photo-induced force microscopy (PiFM) measures photo-induced near-field attractive force to acquire infrared absorption spectrum from nanoscale (~10 nm) region of the sample.





## Agreement between Nanoscale PiFM and Bulk FTIR Spectra

The agreement is excellent on homogeneous samples.



PiFM and FTIR of Nylon













### PVDF PiFM Spectrum Verified via Online Spectral Search Library



PiFM shows excellent agreement with FTIR and is compatible in FTIR spectral search engines.

#### Output from FTIRsearch.com



	Circle of a composite name below to toggie view of a single spectrum, a spito to spectral						
Row	HQI	Compound Name					
1	0.994	POLYVINYLIDENE FLOURIDE	<u>Details</u>				
2	0.993	Poly(vinylidene fluoride)	<u>Details</u>				
3	0.992	CHROMATE(VI); CESIUM; HYDRATE	<u>Details</u>				
4	0.992	Alkylaryl sulfonate	<u>Details</u>				
5	0.992	Alkylaryl sulfonate	Details				
6	0.992	Petroleum sulfonate, Na salt + mineral oil	<u>Details</u>				



#### PCL/5AA – Crystals Topography and PiFM Spectra (Example of PiF-IR Spectra on Inhomogeneous Sample)



#### Exceptional Spatial Resolution in Chemical Mapping









A-A': 21 nm B-B': 6.5 nm



## SIS AlO<sub>x</sub> PS-b-PMMA (41 nm FP Lamellar) PiFM Images



#### Nanoscale Defect Identification via PiF-IR Spectra



## Effect on ODTS Packing Density due to the Nature of Oxide



## Trace Residue Visualization & Identification





#### PiFM Sensitivity - Incomplete Monolayer of Peptoid Molecules



## Complements other Analytical Techniques

	IR PiFM	Raman	FTIR	TOF-SIMS	XPS	TXRF	SEM/EDS	TEM	Auger
Species Detected	M.I.	M.I.	M.I.	M.I.	M.I.	E.I.	E.I.	E.I.	E.I.
Imaging/Ma pping	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Lateral Resolution	~ 5 nm	>0.5 μm	> 10 µm	> 0.2 µm	10 μm – 2 mm	~ 10 mm	1 nm* 0.5μm EDS	0.2 nm	> 10 nm
Depth Probed	20 nm & bulk	> 500 nm	1 µm	1 nm	10 nm	10 nm	1 µm	~ 100 nm	10 nm

\* Imaging

M.I. Molecular information

n E.I. Elemental information



# Summary

- IR PiFM is introduced.
- PiFM measures sample's IR absorption via mechanical force detection.
- It provides for exceptional spatial resolution (~ 5 nm), excellent surface sensitivity (monolayer), and ease-of-use.
- Vista-IR is a turnkey PiFM systems with visible to mid-IR lasers.
- Vista-IR can chemically map and identify organic and inorganic defects and residues via localized IR spectrum (from < 10 nm x 10 nm region by monolayer volume).</li>

