#### HORIBA Scientific Particle Characterization

# Exploiting Fluorescence for Enhanced Nanoparticle Tracking Analysis of CMP Slurries

Jeffrey Bodycomb, Ph.D.



October, 2021

### **Background – SEMI Standards**

Semi standards do not list FCS or NTA on the list of techniques for particle size...

However, they do acknowledge the needs for improved particle metrology as semiconductor technology advances.

#### SEMI C98-1219

### GUIDE FOR CHEMICAL MECHANICAL PLANARIZATION (CMP) PARTICLE SIZE DISTRIBUTION (PSD) MEASUREMENT AND REPORTING USED IN SEMICONDUCTOR MANUFACTURING

#### "Response to Semiconductor Technology Advances

As node sizes continue to decrease, new requirements for CMP PSD are becoming increasingly important.

- As CMP slurries move towards smaller particle sizes more nano-size sensitive PSD metrologies are needed.
- More PSD information is needed at smaller particle sizes (i.e., more particle size bins in the nanometer size range).
- Metrologies capable of differentiating the nature of the measured particles are needed...."

Next gen PSD requires: "Sensitivity in the nm range....Smaller bins in the nm range"



### **The Problem**

What happens when mean particle size approaches the limit of resolution?



Improve Process Control and meet Roadmap Challenges: Improve the resolution of fine particles



# **Nanoparticle Tracking (NTA)**



© 2021 HORIBA, Ltd. All rights reserved.

# Why three colors?





Number, not volume based distribution.

Particle concentration!





### Repeatability

### **Dilute and run PL-7 15 times**

Concentration: 8.48E+13 p/mL CoV (SD/mean): 3.54%

Median Size: 120.93 nm CoV (SD/mean): 1.04%

What is driving repeatability?





Do we follow counting (Poisson Statistics)? Zirconia, run six times as repeats





### **Error bar is standard deviation of 6 repeats**





### Yes

Plot standard deviation as a function of number of particles in each bin. Dashed line is sqrt(num particles) and theoretical.

### This means

- 1) mixing is correct (random sampling)
- 2) uncertainty can be estimated by looking at particle count.



### **Volume normalization**



Effective volume goes down as particle size goes down.

This is in surmountable...but noise (uncertainty due to fewer particles) is not.

### Look to fluorescence to help



Scatter and fluorescence measured orthogonally.

With 3 laser excitation sources, we have many options of fluorophores to selectively identify mixed particles systems. Long pass and short pass filters allow for isolation of fluorescent particles.



### **Detection limits**



in water (n=1.337), laser 445 nm wavelength, polarized, objective NA=0.28, 80°÷100° integration

	Rayleigh Scattering		
$I=I_0\left(rac{1+\cos^2 heta}{2R^2} ight)\left(rac{2\pi}{\lambda} ight)^4\left(rac{n^2-1}{n^2+2} ight)^2\left(rac{d}{2} ight)^6$			
	I	=	scattered intensity
	Theta	=	viewing angle
	R	=	viewing distance
	lamda	=	wavelength of light
	n	=	particle refractive index

particle refractive index = particle diameter =

Rayleigh scattering diminishes quickly ~ d<sup>6</sup>

d

Calculated limit for detecting scatter in silica is ~ 40 nm.

Fluorescent particle emission signal diminishes less rapidly,  $\sim d^2$ 



## **Enhancing detection limits**



combine the scattered signal with fluorescent emission



## **Test the Concept: A Challenging Colloidal Silica**



D2 = 40 nm



With a ~ 40nm mean particle size and a wide range of aggregate sizes, Fuso PL-1 is a good challenge for the resolution of the Viewsizer 3000.



https://fusokk.co.jp/eng/wp/wp-content/uploads/2019/05/5d5d380df5aea451f62d038624cb8cb5.pdf

HORIBA

Scientific

# **Test the Concept: An Efficient Dye**

Cationic dye can interact electrostatically with the negative surface charge colloidal silica

• Very efficient adsorption in the literature

Surface modification of colloidal silica particles using cationic surfactant and the resulting adsorption of dyes *Asad M. Khan* Journal of Molecular Liquids, Volume 274, 15 January 2019

Selective adsorption of organic dyes by porous hydrophilic silica aerogels from aqueous system *Wei Wei* **Water Science & Technology | 78.2 | 2018** 

Adsorptive removal of methylene blue from aqueous solution using coal fly ash-derived mesoporous silica material *Ning Yuan* **Adsorption Science & Technology 2019, Vol. 37(3–4) 333–348** 

Rhodamine 6G (R6G)  $\lambda_{ex}$  = 530 nm,  $\lambda_{em}$  = 565 nm

Based on Prof. Remsen's results and literature supporting cationic dye adsorption on SiO2, and good fluorescence alignment with our green laser We chose a high purity colloidal silica **and Rhodamine 6G** 



# **Test the Concept: Series of Experiments**



Diameter, nm

# **Optimized Conditions**



# **Series of Experiments**



Comparing the calculated limit of scattering detection for the Viewsizer and the enhanced signal of the dyed particles, The improved detection is pronounced at sizes smaller than ~ 65nm.

### **Illustration of the Experiments**



Standard process Fines poorly detected



Dyed and diluted, we can optimize the optical settings to improve detection of fines

Dyed particles increase counts Need dilution to count all particles



# **Summary of Results**

Addition of Rhodamine 6G to Fuso PL-1:

- Consistently increases the counts
- Detection of finer particles upon optimal dilution

CS

More counts->better statistics

Opportunities for follow-up

More efficient dyes

RIBA

- blue and red
- Lot-to-lot (good vs. bad)
- Dyes for other particle types
- Process for complex slurry formulations (interfering additives)





